

# What's new in ProteoIQ version 1.3

The latest version of ProteoIQ – 1.3 - contains the following key enhancements and features listed below.

Main Features
<b>BioWorks .SRF</b> results can now be imported into ProteoIQ
<b>Normalization by selected proteins</b> as internal controls
<b>Combine Mascot, SEQUEST and X!Tandem</b> results into a single ProteoIQ project
<b>Option to export based on protein set view</b>
<b>Filter by database search engine type</b>
<b>Customize resolution</b> on all exported charts

## Downloading and Installing

### System Requirements

Before installing, make sure your computer meets the system requirements. ProteoIQ **MUST** be installed on a computer with **minimum of 2 GB** of physical memory (RAM). It is also recommended that the computer have at least 4 GB of free disk space on the harddrive containing the ProteoIQ installation.

### Downloading ProteoIQ Version 1.3

To access new software versions, log into your account at <https://www.bioinquire.com/signin-register.php>

Please note that the software download links listed in your account page are for 32 bit systems **ONLY**. If you require a 64 bit installation please contact BioInquire.

## Backup ProteoIQ Files

Prior to installing version 1.3 you **MUST** backup the following files:

1. Copy the **db**s folder located at C:\Program Files\ProteoIQ\db to another location. Once the installation is complete the db
s folder should be copied back into the path above.2. Create a backup copy of the file **serial.dat** file located at C:\Program Files\ProteoIQ\

## Removing older versions of ProteoIQ

**Please Note** that you must first **Uninstall** your previous installation of ProteoIQ prior to installation of Version 1.2. Prior to uninstalling it is recommended that you back up the dbs folder located at C:\Program Files\ProteoIQ\db and the serial.dat file located at C:\Program Files\ProteoIQ\. After installing version 1.3 copy the contents of the dbs folder and the serial.dat file back into the ProteoIQ installation directory.

## Installing ProteoIQ Version 1.3

Double click on the ProteoIQ installer and follow the instructions to install ProteoIQ. For **Windows Vista** users, you may need to right click on the ProteoIQ installer and select “Run as administrator”.

## For Technical Assistance:

Email: [support@bioinquire.com](mailto:support@bioinquire.com)

Phone: 706-583-4000

Fax: 706-583-4037