

What's new in ProteoIQ version 2.0

The latest version of ProteoIQ – 2.0 - contains the following key enhancements

Main Features

Biological Annotations Speed Discovery of Biological Significance

- Databases preloaded with GO Terms and KEGG Pathways
- Annotate your own databases with custom annotations

ANOVA Statistic for Expression

- Assess the significance in protein expression changes

Improved and Expanded Normalization

- More control over data normalization including apportionment of shared peptides

Greater Ease of Use

- More Filters
- Check boxes making selection/removal of proteins and peptides easier than ever

More File Formats Supported

- ProteomeDiscoverer (.msf)
- Sequest zipped (.out)
- Generic (.pepXML)

Performance Optimization

- Improved memory management making projects size smaller
- New spectral storage increases speed in project transfers and compression

Downloading and Installing

System Requirements

Before installing, make sure your computer meets the system requirements. ProteoIQ **MUST** be installed on a computer with **minimum of 2 GB** of physical memory (RAM). It is also recommended that the computer have at least 4 GB of free disk space on the hard drive containing the ProteoIQ installation.

Downloading ProteoIQ Version 2.0

To access new software versions, log into your account at:

<https://www.bioinquire.com/signin-register.php>

Installing ProteoIQ Version 2.0

Double click on the ProteoIQ installer and follow the instructions to install ProteoIQ. For **Windows Vista or 7** users, you must install ProteoIQ as administrator (right-click on the installer and/or the ProteoIQ program and select "Run as administrator").

For Technical Assistance:

Email: Kyle.Jones@NuSep.com

Phone: 1-877-592-1060

Fax: 1-877-283-3788